

# HAMA-SURF™ ALCOHOL

## Quaternary Ammonium Compounds for Hard Surface Disinfectant

**READY TO USE**



**For rapid disinfection of all hard surfaces and sensitives in hospital including worktops, hand instruments etc Specially formulated for use in the medical & dental environment.**

### Effective spectrum

Methicillin Resistant Staphylococcus Aureus (MRSA), H1N1 Influenza A Virus, E-Coli, Pseudomonas aeruginosa, Enterococcus hirae, Mycobacterium tuberculosis (TB), Clostridium Difficile (C-Diff) vegetative cell formation (growing cells) of Gram positive organisms, HIV-1, Human Influenza B Virus (HIBV), Staphylococcus Aureus, Hepatitis C Virus (HCV) Aspergillus niger, Candida albicans.

### Ready to use

Spray surfaces and objects to be disinfected from a distance of approx. 30 cm so that they are completely moistened and let the disinfectant take action over a period of 30 seconds. Allow surface spray to evaporate or wipe with a clean tissue or surface wipe.

### Composition

Aqua water, Didecylidimethylammonium Chloride, Ethylenediaminetetracetic Acid, Ethanol, Parfum, Limonene and Non ionic surfactant

### Physical Properties

Light Yellow liquid Odour Lemon, PH 8-9 of solution with a shelf life of 3 years.

### Packaging

1111-132 750 ml bottle  
1111-121 1000 ml bottle



Lemon Aroma

